


SigMate Installation Guide

V 0.1 (Beta)


For any suggestion, question, request to adapt the code
to your analysis requirement, please contact:

muftimahmud@gmail.com


- Extract the content of the SigMate.zip file.




- Put it in the Matlab's working directory.




- Browse to the SigMate directory in the Matlab console's 'Current Directory'.



- Type 'FirstGUI' (without quotes) in the console's command window.



- Typing 'FirstGUI' will open the main GUI of SigMate from where using the menubar the other features can be accessed.



- Important Notes:
 - when the initial GUI is opened, the SigMate folder and its first level subfolders are added in the path variable, and closing this GUI will remove the folders from the path.
 - each feature, when accessed from the menubar will put its subfolders in the path variable and closing it will remove them from the path.